AT52323

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ENG

INTRODUCTION

Thank you for purchasing our product. Before using, please read this user manual carefully. Incorrect use may lead to health detriment, injury or accidents. Therefore, in order to use the product safely and successfully, please read this manual carefully.

DESTINY

The electric wheelchair is a new, original product that will provide excellent travel comfort to people with mobility difficulties. An innovative and intelligent wheelchair will help people with disabilities in their daily activities.

CHARACTERISTICS

The ergonomic shape and the possibility of personalizing the stroller provide high comfort of use.

REMARKS

This electric wheelchair has been tested according to ISO 7176-9 Climatic Testing.

The study guarantees that people using the electric wheelchair will be able to find shelter in the event of precipitation.

IT IS NOT ALLOWED to use the electric wheelchair during an electrical storm;

IT IS NOT ALLOWED to expose the electric wheelchair to rain;

It is FORBIDDEN to store the electric wheelchair in a damp place for a long time;

It is NOT ALLOWED to spray, wash or clean the electric wheelchair with an automatic car washer.

Direct exposure to rain, smoke or humid air may cause electrical and mechanical failure of the electric wheelchair as well as lead to the formation of rust.

TERMS OF USE

- Ambient temperature: -10 ° C ~ 50 ° C; Relative Humidity: 20% ~ 90%;
- Atmospheric pressure: 80kpa ~ 1.23kpa; No strong electromagnetic interference

USE

1. Navigating on slopes

This electric wheelchair has successfully passed tests on inclines up to a maximum angle of 10 ° with a load of 120 kg. NEVER drive on slopes greater than 10 ° Restrict the speed of the wheelchair when backing up a slope. NEVER back up on slopes greater than 10 °. Also, take extra care when backing up.

Slope performance and range depend on battery and usage conditions:

For example: The weight of the person using the wheelchair; Surface (e.g., grass or gravel); Tilt; Battery life; Extreme temperatures; Application and weight of additional equipment etc.

2. Brake

When braking, maintain the correct position, grasp the handrail and release the control handle. The electric wheelchair will stop in a few seconds.

Note: The electromagnetic brake does not activate immediately, but only after 1/2 revolution of the wheel.

3. Emergency braking

In the event of an emergency or unexpected movement, the wheelchair can be stopped using the power switch (this method of braking, however, should not be used under normal conditions of use).

Note: Frequent use of the emergency brake may damage the motor, which in turn causes the electric wheelchair to malfunction.

4. Using a mobile phone

It is FORBIDDEN to use cell phones or wireless devices while operating the electric wheelchair.

The use of cell phones or wireless devices generates strong electromagnetic fields which interfere with the electronic system of the electric wheelchair. If you need to use cell phones or wireless devices, the electric wheelchair must be stopped and turned off beforehand.

WARNINGS

IT IS NOT ALLOWED to bend down to pick up objects while in the electric wheelchair; such action may overturn the wheelchair and endanger the health of the user.

IT IS FORBIDDEN to use the electric wheelchair on the escalator as this may cause a risk of injury;

IT IS FORBIDDEN to move the wheelchair in the lane of vehicles and off-road;

NEVER move on surfaces covered with water, oil, ice, etc. Driving on slippery surfaces may cause you to lose control of the wheelchair;

IT IS FORBIDDEN to try to drive over obstacles and ditches exceeding 45 mm, as it may cause bodily injuries;

NEVER make sudden maneuvers or reverse the wheelchair at excessive speed;

IT IS FORBIDDEN to use the trolley to transport people;

NEVER tow any other electric wheelchair;

It is FORBIDDEN to use the electric wheelchair without checking the correct connection and fastening of the fasteners. Always pay attention to your surroundings when using your wheelchair;

DO NOT make unauthorized repairs to parts, accessories or adapters of the electric wheelchair;

You MUST NOT get on the wheelchair without applying the brake;

IT IS FORBIDDEN to connect any other medical devices (such as life support devices, ventilator, etc.) to the electronic system of the electric wheelchair;

IT IS FORBIDDEN to use the electric wheelchair if you have a cognitive impairment (e.g. after drinking alcohol, taking drugs or medications);

IT IS FORBIDDEN to use the electric wheelchair in the event of a sudden loss of consciousness;

IT IS FORBIDDEN to use the electric wheelchair if your eyesight is impaired; The manufacturer declines any responsibility for accidents resulting in personal injury caused by the violation of local or international laws or regulations; The stroller may only be used after properly seated on the chair;

IT MUST NOT rest your feet on anything other than the footrests;

IT IS NOT ALLOWED to stand on the seat of the electric wheelchair;

It is FORBIDDEN to attempt to ride up or down slopes greater than 10 °;

IT IS NOT ALLOWED to back up on slopes greater than 10 °; When driving on slopes, release the control lever. In this situation, the electric wheelchair may reverse. When moving forward or backward, the electric wheelchair will start to brake after driving a distance of approximately 30 cm; Before using the electric wheelchair, check the connection of all electrical parts for safety and reliability;

It is NOT ALLOWED to disconnect, cut or modify the wiring harnesses installed on the electric wheelchair.

The use of unauthorized batteries such as non-recyclable lead acid batteries etc. is NOT ALLOWED.

Before installing, read the information on the batteries and charger. Violation of the above requirements may create unexpected threats.

DESCRIPTION

- Energy-saving and efficient DC brush motor is the main driving part of the electric wheelchair. The motor is characterized by high torque, high efficiency and energy saving.
- The advantages of the DC brush motor are high efficiency, long life and low energy consumption. The motor noise can be controlled within 65dB.
- The use of a lithium battery as a green energy source in an electric wheelchair has advantages such as high energy density, long life, energy saving and environmental protection.
- Dedicated brush driver enables intelligent motor control. The control of the throttle allows you to drive forward, turn and reverse while limiting the speed.
- Patented technologies are used in the frame structure. The folding mechanism based on four crossbars allows for quick and effective folding and unfolding. The electric wheelchair is equipped with an electric or manual folding function, and has the advantage of safety, reliability and ease of use.
- The foldable structure facilitates the transport of the stroller and allows it to be folded to a small size (to be placed, for example, in the trunk or for transport by air).
- The gyro position detection function allows you to get a voice alert for unevenness and tilt. Reducing the speed of the electric wheelchair to gear 1 will trigger an alarm in the app
- If you do not perform any operations (no operation of buttons, levers or APP) within 30 minutes, the electric wheelchair will automatically shut down when it is in energy saving mode
- To set a constant cruising speed, hold down the "speed +" button for a longer time, then push the lever forward to activate the cruise control
- Wireless control. The wireless remote control is equipped with a lever that allows you to move forward, reverse, turn left and right. In addition, a dedicated button allows you to automatically fold and unfold the electric wheelchair.
- Beginner work mode

GENERAL INFORMATION

SUSPENSION

The stroller is equipped with an independent suspension system. The front wheels have two sets of springs that absorb unevenness, stabilize the stroller during the ascent and provide high comfort.

TECHNICAL PARAMETERS

Product model: AT52323

Overall dimensions: 1075mm long, 628mm wide, 930mm high

Storage dimensions: 895mm long, 628mm wide, 395mm high

Handle height: 280 ± 5mm

Seat height: 52 cm

Backrest height: 45 cm

Seat width: 42 cm

Maximum user weight: 120kg

Folding Mode: Electric

Maximum speed: 7 km / h

Braking performance on slopes: 12 °

Maximum height of the barrier that can be overcome: 45mm

The maximum width of the step that can be overcome: 100mm

Gradeability: 10 °

Own weight: 42.1 ± 2 kg

Battery Quantity: 1

Battery weight: 3.4 ± 0.15 kg

Solid front wheels - 10 inch multi-directional;

Rear wheels - 10 inches, PU tire (2.4MPa-2.7MPa)

Motor parameters: DC motor rated voltage 24V Rated power 200W * 2

Lithium battery parameters: DC24V20AH

Charger parameters: The input voltage is 29.4V and the output current is 4A

Static stability: 20 ° for the longitudinal direction and 15 ° for the lateral direction

Dynamic stability: 10 ° for both climbs and slopes

Braking distance on a flat surface: 100 cm

The maximum output current of the control system: 30A

Approximate range on a single charge: 25 km

WATERPROOF LED REAR LIGHTING

Suspension: independent suspension

Ground clearance height: 60mm

Minimum turning radius: 730 mm


Waterproof level: IPX4

Pedal angle adjustment: adjustable

STORAGE

- Store the packed electric wheelchair at -40 ° C-60 ° C
- The humidity level cannot exceed 95%
- Atmospheric pressure - 56-146kPa
- The electric wheelchair should be stored in a room not exposed to harmful gases that may cause corrosion, provide adequate ventilation. Avoid excessive loading and stacking.
- When stacking is required, it must be stored in a maximum of three layers.
- Do not store the product for a long time in places near high temperature sources or exposed to direct sunlight.
- When storing the trolley for an extended period, turn off the power switch after it is fully charged.
- If the product is stored for more than one month, fully charge it, then continue storing and using it.

DESCRIPTION OF COMPONENTS



USE

- When using the electric wheelchair for the first time, you must move on flat and solid ground, free from obstacles and bystanders.
- The maximum load of the electric wheelchair is 120 kg, it is forbidden to drive under excessive load; the trolley may only be used by 1 person;
- Before use, make sure that the electric wheelchair is fully unfolded;
- You must also fasten your seat belt;
- Correctly sit on the electric wheelchair, slightly push the control lever forwards. The electric wheelchair will move forward and backward and turn left or right; when released, the electric wheelchair will stop. Practice these basic functions until you are proficient in the operation;
- When using the wheelchair, remove any nearby obstacles;
- NEVER drive the wheelchair towards bystanders;
- Reduce the speed of the truck when complex maneuvers (eg driving through doors or steering) must be performed. In severe conditions, the speed of the electric wheelchair should be limited to the lowest level;
- Before using the trolley, make sure that the device is working properly;
- The electric wheelchair switches itself off automatically after 30 minutes of inactivity.

UNFOLDING AND FOLDING

While folding or unfolding, deflecting the handle will stop folding and unfolding. Stopped folding or unfolding process can be continued as described in steps 1 and 2


Press the "horn button" and "speed -" simultaneously on the electric wheelchair for 5 seconds

Press the "horn button" and "speed +" on the electric wheelchair simultaneously for 3 seconds.

The unfolding and folding of the stroller will be completed automatically.

CONTROL

1. Electric battery charge indicator
2. Speed
3. Speed gearbox
4. Speed acceleration
5. Sound signal
6. Reduction in speed
7. Power / Lockout
8. Lockout status indicator light
9. Button for starting the work mode for a beginner
10. SOS button
11. Shifter
12. Displaying the battery% value
13. Bluetooth status display
14. Error code display

Shifter

The direction of travel and speed of the electric wheelchair are controlled by pressing the control lever;

Power switch / lock key

2 seconds long press turn on the power, long press 2 seconds turn off the power;

Pressing the power button will turn off the wheelchair's electric lock indicator and the wheelchair will unlock

Press the power button, the electric wheelchair lock indicator will be displayed and the electric wheelchair will be locked

"Speed +" button

Pressing the button increases the speed;

* To set a constant cruising speed, hold down the "speed +" button for a longer time, then push the lever forward to activate the cruise control

"Speed-" button

Pressing the button reduces the speed;

Beep button

A short press of the beep button will trigger the beep. Sound intensity meets national noise standards;

Electric trolley lock status indicator light

Turning on the power automatically activates the lock. Re-activation of the lock is possible at any time while using the device;

When the lock is engaged, the wheelchair will not respond to movement of the joystick. When the lock is released, deflecting the control stick will cause the wheelchair to move.

When using the application, information about the lock status will be displayed.

SOS emergency button

Holding the button down for 3 seconds will trigger an alarm signal

Start-up button for a beginner

The driving parameters of the electric wheelchair will be automatically adjusted to the speed and parameters adapted to the learners.


Switching the truck to manual mode

Press the "+" button and the "-" button simultaneously to put the truck into manual drive mode. Move the control lever to exit manual drive mode.

BATTERY

Removing the battery

- 1 - Open the battery lock located on the housing of the trolley
2. Pull out the battery

MAINTENANCE AND CLEANING

- We recommend regular (weekly or monthly - depending on frequency of use) inspection of fasteners such as bolts and nuts for possible loosening, possible corrosion etc. ;
- Regularly check the moving parts of the joint for binding, falling, etc. ;
- Regularly check the tires for signs of wear, etc;
- When using a wheelchair in rain or humid weather, wipe the device dry as soon as possible to prevent rust formation;
- Avoid exposing the electric wheelchair to excessive sunlight in high temperature conditions.
- Cleaning the electric wheelchair: Wipe the surface of the wheelchair with a clean, soft cloth.

TROUBLESHOOTING

When calling the service, have the fault code ready.

The product is equipped with an automatic diagnosis system, thanks to which a fault code is displayed

The fault code	The nature of the fault	
E1	Button malfunction	Emergency signal
E2	Lever malfunction	Emergency signal
E3	GPRS failure	Emergency signal
E4	Battery failure	Emergency signal
E5	Relay error	Emergency signal
E6	Accessory fault	Emergency signal
E7	Bluetooth malfunction	Emergency signal
E8	Engine error	Emergency signal
E9	Drive error	Emergency signal
EA	FLASH failure	Emergency signal

Warning text	Workflow
<i>Enter novice mode</i>	Press the start button for the beginner mode
<i>Quit novice mode</i>	Press the start button for the beginner operation again
<i>Alarm</i>	Long press of the SOS button
<i>Lock</i>	By default, your device starts in Locked Mode and then goes from Unlocked to Locked Mode
<i>Unlock</i>	Unlocking the lock
<i>Wireless remote mode</i>	Operation via the control button on the remote control
<i>Quit remote control mode</i>	Operation via the control button on the remote control
<i>Fold</i>	Simultaneous long press of the beep button and the SPD button
<i>Expand</i>	Simultaneous long press of the beep button and the SPD + button
<i>Upgrade mode</i>	After the application detects a new version, press the button to install
<i>Upgrade finished</i>	It will be launched after updating a single module
<i>Update</i>	The device is updating
<i>Attention, reserve</i>	Low charge warning signal
<i>Power on sound</i>	Power on, the power on signal will sound after 3 seconds
<i>Charging</i>	<i>Charging</i>
<i>Charging completed</i>	After charging is complete
<i>Band-type brake open</i>	The wheelchair brake is deactivated
<i>Power off sound</i>	Power off, there will be a power off signal after 3 seconds
<i>Equipment bumpiness</i>	Moving on uneven surfaces
<i>Equipment slant</i>	The machine will be on an incline of more than 5 °
<i>Sharp bend</i>	The signal will be triggered when there is a sharp turn
<i>Upgrade failed</i>	Update error
<i>Version V1.0</i>	Voice prompts version
<i>Lower power</i>	Battery life is less than 20%, it will announce every minute
<i>Stop Electric folding wheelchair</i>	<i>Stop Electric wheelchair folding</i>
<i>Stop Electric unfolding wheelchair</i>	<i>Stop Electric unfolding wheelchair</i>
<i>Pushing mode</i>	Press the "+" and "-" keys simultaneously
<i>Quit pushing mode</i>	Pressing the lever exits the manual drive mode

CAUTION:

When a product-related "serious incident" occurs, which directly or indirectly led, could or may have led to any of the following events:

- a) death of a patient, user or other person or
- b) temporary or permanent deterioration in the health of a patient, user or other person, or
- c) a serious threat to public health

this "serious incident" must be reported to the manufacturer and to the competent authority of the Member State in which the user or patient is resident. In the case of Poland, the competent authority is the Office for Registration of Medicinal Products, Medical Devices and Biocidal Products.

CAUTION:

In the event of pain, allergic reactions or other disturbing, unclear for the user symptoms related to the use of a medical device, please consult a healthcare professional.

ATTENTION: it is forbidden to use the product other than in accordance with its intended purpose

ATTENTION: The manufacturer is not responsible for damage caused by neglect of maintenance, inadequate servicing or resulting from failure to comply with the recommendations contained in this manual.

NOTE: There may be a risk of tipping over if the product is used improperly. Please follow the boarding / dismounting / moving instructions.

ATTENTION: During use and handling of the product, as well as during its assembly and adjustment of mechanisms, there may be a risk of entrapment and / or pinching of body parts of the user / accompanying person in the openings / gaps between the elements. These steps must be performed with extreme caution. After completing the adjustment, stabilize the position by tightening the nuts / bolts carefully

This symbol indicates the user's maximum weight

METHOD OF DISPOSAL OF THE PRODUCT AFTER REMOVING FROM SERVICE

After the device is withdrawn from use, the medical device can be disposed of as normal household waste, except for electrical products - proceed in the manner appropriate for the disposal of electrical and electronic equipment.

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Model:
Serial number:

Point of sale stamp

and a legible signature of the Seller:

1. AN TAR Sp. J., 03-068 Warsaw, ul. Zawiślańska 43, hereby grants a 12-month warranty on the rollator, starting from the date of purchase by the Buyer. Products financed by the National Health Fund are covered with the extended warranty period, which covers half of the useful life specified in the Regulation on Guaranteed Services, provided that the useful life is specified in the Regulation in months or years.
2. During the warranty period AN TAR undertakes to carry out the repairs necessary to restore proper functioning of the product, within 14 days from the date of receipt of the product by the service centre, together with the proof of purchase and the warranty card.
3. AN TAR reserves the right to replace the product if the repair is deemed uneconomic. Repair or replacement of faulty products does not extend the warranty period.
4. If the product is found to be malfunctioning, send it immediately to AN TAR (as specified in point 1) or contact the point of sale.
5. Before use, read the instructions carefully and follow them. If the product is used improperly and contrary to recommendations, the warranty will be irrevocably lost.
6. The warranty does not cover damage caused by external mechanical force, load exceeding the maximum permissible weight of the user as described in the manual, as well as contact of the product with alcohol, greasy substances or petrol. The warranty does not cover parts subject to normal wear and tear during use, such as the bearing, rubber parts and others.
7. The product delivered to the service in connection with the complaint cannot cause an epidemic.
8. THIS WARRANTY IS ONLY VALID TOGETHER WITH THE PROOF OF PURCHASE (BILL, RECEIPT, INVOICE). IN CASE OF A COMPLAINT, A PROOF OF PURCHASE MUST BE ATTACHED. IN ORDER TO EXTEND WARRANTY PERIOD FOR PRODUCTS FINANCED BY THE NATIONAL HEALTH FUND, ATTACH A COPY OF THE ORDER FOR THE SUPPLY OF MEDICAL DEVICES BEING ORTHOPAEDIC ITEMS AND AIDS.
9. Complaints submitted without the proof of purchase and a warranty card with the product serial number will not be considered.

 AN TAR Sp. J.
03-068 Warsaw, ul. Zawiślańska 43
Tel. 22 518 36 00, Fax 22 518 36 30
www.antar.net

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